



Fiber Optic Cannulae

- 80% Efficiency cut-off
- Black Ceramic ferrule ensures low autofluorescence
- Cannulae optimized for fiber photometry use
- Ferrule sizes available:
1.25 mm or 2.5mm
- Various Numerical Aperture (NA) available: 0.22, 0.37, 0.39, 0.48, 0.50
- Standard fiber length = 10 mm
Custom lengths (2mm-40mm) can be made to order
- Available in standard core size of 200 or 400 um.
- Bulk discounts available for fiber orders of 100+ fiber optic cannulae

